

No Corners, No Loss.

Trust Stericup® and Steriflip® filters with Millipore Express® PLUS membrane to ensure fast and complete recovery, with no loss of media in the filtration process.

SPECIAL OFFER

Buy 3 cases in any combination of Stericup®, Steritop® & Steriflip® filters

GET ONE CASE FREE!*


Use promotion code ZBIOS651

To learn more about sterile filtration products, visit:
www.millipore.com/cellculture

Offer expires December 31, 2013.



Choose Stericup® filter units for fast, complete recovery:



Stericup®, Steritop® and Steriflip® Filters: The Industry Standard for Sterile Filtration Quality and Efficiency

Stericup® Filter Units – combine a filter unit with a receiver flask and cap for processing and storage.

Description	Membrane	Pore Size (µm)	Funnel Capacity (mL)	Receiver Bottle (mL)	Qty/Pk	Catalogue No.
Stericup®-VP Filter Units	Millipore Express® (PES)	0.1	250	250	12	SCVPU02RE
			1000	1000	12	SCVPU11RE
Our featured 250 mL and 1000 mL Stericup® units provide fast flow and low-protein-binding 0.1 µm PES membrane for sterilization and mycoplasma removal!						
Stericup®-GP Filter Units STEM CELL TESTED	Millipore Express® PLUS (PES)	0.22	150	150	12	SCGPU01RE
			250	250	12	SCGPU02RE
			500	500	12	SCGPU05RE
			500	1000	12	SCGPU10RE
			1000	1000	12	SCGPU11RE
Stericup®-GV Filter Units	Durapore® (PVDF)	0.22	150	150	12	SCGVU01RE
			250	250	12	SCGVU02RE
			500	500	12	SCGVU05RE
			500	1000	12	SCGVU10RE
			1000	1000	12	SCGVU11RE
Stericup®-HV Filter Units	Durapore® (PVDF)	0.45	150	150	12	SCHVU01RE
			250	250	12	SCHVU02RE
			500	500	12	SCHVU05RE
			1000	1000	12	SCHVU11RE

Steritop® Bottle Top Filter Units – can be used with bottles having a 33 mm or 45 mm opening.

Description	Membrane	Pore Size (µm)	Funnel Capacity (mL)	Receiver Bottle (mL)	Qty/Pk	Catalogue No.
Steritop®-GP Filter Units STEM CELL TESTED	Millipore Express® PLUS (PES)	0.22	150	33	12	SCGPS01RE
			150	45	12	SCGPT01RE
			250	33	12	SCGPS02RE
			250	45	12	SCGPT02RE
			500	33	12	SCGPS05RE
			500	45	12	SCGPT05RE
			1000	45	12	SCGPT10RE
Steritop®-GV Filter Units	Durapore® (PVDF)	0.22	500	45	12	SCGVT05RE

Steriflip® Filter Units – completely eliminate liquid handling for sterile samples you can trust.

Description	Catalogue No.
0.22 µm Millipore Express® (PES) Steriflip®-GP Filter unit	SCGP00525
0.45 µm Durapore® (PVDF) Steriflip®-HV Filter unit	SE1M003M00

*Order any three packs of Stericup®, Steritop® or Steriflip® filters and get **ONE FREE** filtration product of equal or lesser value. To receive your free product, use promotion code **ZBIOS651**. Offer valid until December 31, 2013. **Limit: one order per customer.** This promotion (**ZBIOS651**) is for one-time-use only and cannot be duplicated or transferred.

Get Connected!

Join Merck Millipore Bioscience on your favorite social media outlet for the latest updates, news, products, innovations, and contests!

 facebook.com/MerckMilliporeBioscience

 twitter.com/Merck4Bio



www.merckmillipore.com

Merck Millipore and the M mark are trademarks and Stericup, Millipore Express, Durapore, Steritop, and Steriflip are registered trademarks of Merck KGaA, Darmstadt, Germany. Lit No. PR5338ENEU 01/2013 BS GEN-13-07732 ©2013 EMD Millipore Corporation, Billerica, MA, USA.

To Place an Order or Receive Technical Assistance

In Europe, please call Customer Service:

France: 0825 045 645

Germany: 01805 045 645

Italy: 848 845 645

Spain: 901 516 645 Option 1

Switzerland: 0848 645 645

United Kingdom: 0870 900 4645

For other countries across Europe, please call: +44 (0) 115 943 0840

Or visit: www.merckmillipore.com/offices

For Technical Service visit:

www.merckmillipore.com/techservice